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1. Zhu K, Guo C, Lai H, Yang W, Wang C. Novel hyperbranched polyamidoamine nanoparticle based gene delivery: transfection, cytotoxicity and in vitro evaluation. *Int J Pharm.* 2012 Feb 28;423(2):378-83.
2. Ou J, Pan F, Geng P, Wei X, Xie G, Deng J, Pang X, Liang H. Silencing fibronectin extra domain A enhances radiosensitivity in nasopharyngeal carcinoma involving an FAK/Akt/JNK pathway. *Int J Radiat Oncol Biol Phys.* 2012 Mar 15;82(4):e685-91.
3. Shi S, Deng YZ, Zhao JS, Ji XD, Shi J, Feng YX, Li G, Li JJ, Zhu D, Koeffler HP, Zhao Y, Xie D. RACK1 promotes non-small-cell lung cancer tumorigenicity through activating sonic hedgehog signaling pathway. *J Biol Chem.* 2012 Mar 9;287(11):7845-58.
4. Zhu K, Lai H, Guo C, Xu D, Wang C. Novel vascular endothelial growth factor gene delivery system-manipulated mesenchymal stem cells repair infarcted myocardium. *Exp Biol Med (Maywood).* 2012 Jun 1;237(6):678-87.
5. Ou J, Luan W, Deng J, Sa R, Liang H. α V Integrin Induces Multicellular Radioresistance in Human Nasopharyngeal Carcinoma via Activating SAPK/JNK Pathway. *PLoS One.* 2012;7(6):e38737.
6. Zhu K, Guo C, Xia Y, Lai H, Yang W, Wang Y, Song D, Wang C. Transplantation of novel vascular endothelial growth factor gene delivery system manipulated skeletal myoblasts promote myocardial repair. *Int J Cardiol.* 2013 Oct 3;168(3):2622-31.
7. Li H, Qin J, Jin G, Zou L, Shi J, Han X, Cheng X, Zhang X. Overexpression of Lhx8 inhibits cell proliferation and induces cell cycle arrest in PC12 cell line. *In Vitro Cell Dev Biol Anim.* 2015 Apr;51(4):329-35.
8. Xu R, Bi C, Song J, Wang L, Ge C, Liu X, Zhang M. Upregulation of miR-142-5p in atherosclerotic plaques and regulation of oxidized low-density lipoprotein-induced apoptosis in macrophages. *Mol Med Rep.* 2015 May;11(5):3229-34.
9. Peilin Huang, Baozhen Zeng, Zhuoxian Mai, Juntao Deng, Yueping Fang, Wenhua Huang, Hongwu Zhang, Jinying Yuan, Yen Wei and Wuyi Zhou. Novel drug delivery nanosystems based on out-inside bifunctionalized mesoporous silica yolk-shell magnetic nanostars used as nanocarriers for curcumin. *J Mater. Chem. B.* 2016;4(1):46-56

Version 2018.06.01