碧云天生物技术/Beyotime Biotechnology
订货热线: 400-168-3301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Guo M, Huang T, Cui Y, Pan B, Shen A, Sun Y, Yi Y, Wang Y, Xiao G, Sun G. PrPC interacts with tetraspanin-7 through bovine PrP154-182 containing alpha-helix 1. *Biochem Biophys Res Commun*. 2008 Jan 4;365(1):154-7.
2. Xu TF, Wang XL, Yang JZ, Hu XY, Wu WF, Guo L, Kang LD, Zhang LY. Overexpression of Zip-2 mRNA in the leukocytes of asthmatic infants. *Pediatr Pulmonol*. 2009 Aug;44(8):763-7.
3. Chen H, Wang W, Mo Y, Ma Y, Ouyang N, Li R, Mai M, He Y, Bodombossou-Djobo MM, Yang D. Women with high telomerase activity in luteinised granulosa cells have a higher pregnancy rate during in vitro fertilisation treatment. *J Assist Reprod Genet*. 2011 Sep;28(9):797-807.
4. Zhao W, Zhang L, Yin Z, Su W, Ren G, Zhou C, You J, Fan J, Wang X. Activated hepatic stellate cells promote hepatocellular carcinoma development in immunocompetent mice. *Int J Cancer*. 2011 Dec 1;129(11):2651-61.
5. Zhou CS, Yin ZY, Zhao WX, Zhang L, You JY, Su WX, Fan J, Wang XM. Bear bile inhibits the immunosuppression activity of hepatic stellate cells in vivo. *Hepatogastroenterology*. 2012 Jul;59(117):1529-36.
6. Hu Y, Lu W, Shen Y, Xu Y, Yuan Z, Zhang C, Wu J, Ni Y, Liu S, Cao J. Immune changes of *Schistosoma japonicum* infections in various rodent disease models. *Exp Parasitol*. 2012 Jun;131(2):180-9.
7. Yao SQ, Li ZZ, Huang QY, Li F, Wang ZW, Augusto E, He JC, Wang XT, Chen JF, Zheng RY. Genetic inactivation of the adenosine A(2A) receptor exacerbates brain damage in mice with experimental autoimmune encephalomyelitis. *J Neurochem*. 2012 Oct;123(1):100-12.
8. Zhao W, Ren G, Zhang L, Zhang Z, Liu J, Kuang P, Yin Z, Wang X. Efficacy of mesenchymal stem cells derived from human adipose tissue in inhibition of hepatocellular carcinoma cells in vitro. *Cancer Biother Radiopharm*. 2012 Nov;27(9):606-13.
9. Xue Z, Gao J, Zhang Z, Yu W, Wang H, Kou X. Antihyperlipidemic and antitumor effects of chickpea albumin hydrolysate. *Plant Foods Hum Nutr*. 2012 Dec;67(4):393-400.
10. Liu J, Zhang PS, Yu Q, Liu L, Yang Y, Qiu HB. Kinetic and distinct distribution of conventional dendritic cells in the early phase of lipopolysaccharide-induced acute lung injury. *Mol Biol Rep*. 2012 Dec;39(12):10421-31.
11. Kuang P, Zhao W, Su W, Zhang Z, Zhang L, Liu J, Ren G, Yin Z, Wang X. 18β-glycyrrhetic acid inhibits hepatocellular carcinoma development by reversing hepatic stellate cell-mediated immunosuppression in mice. *Int J Cancer*. 2013 Apr 15;132(8):1831-41.
12. Lun Y, Xia H, Zhang Q, Yu C, Chen N, Li X, Liu S, Lei L. Anti-inflammatory and immunosuppressive activities of 1,3-dicyclopentyl-1,2,3,6-tetrahydropyrimidine-4,5-dicarboxylic acid diethyl ester (ZL-5015). *Int Immunopharmacol*. 2013 Oct;17(2):168-77.
13. Ding F, Yu L, Wang M, Xu S, Xia Q, Fu G. O-GlcNAcylation involvement in high glucose-induced cardiac hypertrophy via ERK1/2 and cyclin D2. *Amino Acids*. 2013 Aug;45(2):339-49.
14. Li S, Li Y, Ning H, Na L, Niu Y, Wang M, Feng R, Liu L, Guo F, Hou S, Chu X, Wang Y, Zhang Y, Zhang H, Huang L, Bi M, Huang Y, Hao L, Zhao Y, Wang C, Wang Y, He Y, Sun C. Calcium supplementation increases circulating cholesterol by reducing its catabolism via GPER and TRPC1-dependent pathway in estrogen deficient women. *Int J Cardiol*. 2013 Oct 3;168(3):2548-60.
15. Zhang Y, Wang X, Wang Y, Liu Y, Xia M. Supplementation of Cyanidin-3-O-β-Glucoside Promotes Endothelial Repair and Prevents Enhanced Atherosclerosis in Diabetic Apolipoprotein E-Deficient Mice. *J Nutr*. 2013 Aug;143(8):1248-53.
16. Li D, Dong H, Li S, Munir M, Chen J, Luo Y, Sun Y, Liu L, Qiu HJ. Hemoglobin subunit beta interacts with the capsid protein and antagonizes the growth of classical swine fever virus. *J Virol*. 2013 May;87(10):5707-17.
17. Fu Y, Zhang Y, Chang X, Zhang Y, Ma S, Sui J, Yin L, Pu Y, Liang G. Systemic immune effects of titanium dioxide nanoparticles after repeated intratracheal instillation in rat. *Int J Mol Sci*. 2014 Apr 22;15(4):6961-73.
18. Tian X, Li E, Yang F, Peng Y, Zhu J, He F, Chen X. In vivo immunotoxicity of SiO₂@(Y_{0.5}Gd_{0.45}Eu_{0.05})₂O₃ as dual-modality nanoprobes. *Int J Mol Sci*. 2014 Aug 7;15(8):13649-62.
19. Chen YC, Wen S, Shang SA, Cui Y, Luo B, Deng GJ. Magnetic resonance and near-infrared imaging using a novel dual-modality nano-probe for dendritic cell tracking in vivo. *Cytotherapy*. 2014 May;16(5):699-710.
20. Ding Y, Liang Y, Deng B, Qiao A, Wu K, Xiao W, Gong W. Induction of TGF-β and IL-10 production in dendritic cells using astilbin to inhibit dextran sulfate sodium-induced colitis. *Biochem Biophys Res Commun*. 2014 Apr 4;446(2):529-34.
21. Zhang M, He J, Hou J, Wu J, Sun M, Cui J, Tian J, Jiang M, Yu B. The immunosuppressant Protosappanin A diminished recipient T cell migration into allograft via inhibition of IP-10 in rat heart transplant. *PLoS One*. 2014 May 5;9(5):e96138.
22. Zhang Q, Xie D, Wang S, You C, Monroig O, Tocher DR, Li Y. miR-17 is involved in the regulation of LC-PUFA biosynthesis in vertebrates: effects on liver expression of a fatty acyl desaturase in the marine teleost *Siganus canaliculatus*. *Biochim Biophys Acta*. 2014 Jul;1841(7):934-43.
23. Zhang T, Yu FL, Yang WX, Ruan MM, Yue ZY, Liu Y, Liu TT, Zhou P, Xia H, Guan JC. Staphylococcal enterotoxin B administration during pregnancy imprints the increased CD4: CD8 T-cell ratio in the peripheral blood from neonatal to adult offspring rats. *J Med Microbiol*. 2015 Jan;64(Pt 1):1-6.
24. He H, Chen K, Wang F, Zhao L, Wan X, Wang L, Mo Z. miR-204-5p promotes the adipogenic differentiation of human adipose-derived mesenchymal stem cells by modulating DVL3 expression and suppressing Wnt/β-catenin signaling. *Int J Mol Med*. 2015 Jun;35(6):1587-95.

25. Zhang C, Wang Q, Xi X, Jiao J, Xu W, Huang J, Lai Z. High serum miR-183 level is associated with the bioactivity of macrophage derived from tuberculosis patients *Int J Clin Exp Pathol*. 2015 Jan 1;8(1):655-9.
26. Lu Y, Liu J, Liu Y, Qin Y, Luo Q, Wang Q, Duan H. TLR4 plays a crucial role in MSC-induced inhibition of NK cell function. *Biochem Biophys Res Commun*. 2015 Aug 21;464(2):541-7.
27. Yang F, Ruan YC, Yang YJ, Wang K, Liang SS, Han YB, Teng XM, Yang JZ. Follicular hyperandrogenism downregulates aromatase in luteinized granulosa cells in polycystic ovary syndrome women. *Reproduction*. 2015 Oct;150(4):289-96.
28. Zhang Q, Di W, Dong Y, Lu G, Yu J, Li J, Li P. High serum miR-183 level is associated with poor responsiveness of renal cancer to natural killer cells. *Tumour Biol*. 2015 Dec;36(12):9245-9.
29. Zhang C, Xi X, Wang Q, Jiao J, Zhang L, Zhao H, Lai Z. The association between serum miR-155 and natural killer cells from tuberculosis patients. *Int J Clin Exp Med*. 2015 Jun 15;8(6):9168-72.
30. Yin L, Zhu Y, Yang J, Ni Y, Zhou Z, Chen Y, Wen L. Adipose tissue-derived mesenchymal stem cells differentiated into hepatocyte-like cells in vivo and in vitro. *Mol Med Rep*. 2015 Mar;11(3):1722-32.
31. Yang J, Ren F, Zhang H, Jiang L, Hao Y, Luo X. Induction of Regulatory Dendritic Cells by *Lactobacillus paracasei* L9 Prevents Allergic Sensitization to Bovine β -Lactoglobulin in Mice. *J Microbiol Biotechnol*. 2015 Oct 28;25(10):1687-96.
32. Jiang S, Xia M, Yang J, Shao J, Liao X, Zhu J, Jiang H. Novel insights into a treatment for aplastic anemia based on the advanced proliferation of bone marrow derived mesenchymal stem cells induced by fibroblast growth factor 1. *Mol Med Rep*. 2015 Dec;12(6):7877-82.
33. Wang YZ, Lv H, Hao YL, Zhang HQ, Li L, Cai GM, Hu M, Jia CX, Feng XG, Kong QX. Suppressive oligodeoxynucleotides induced tolerogenic plasmacytoid dendritic cells and ameliorated the experimental autoimmune neuritis. *Immunol Lett*. 2015 Jul;166(1):13-8.
34. Liu ZJ, Chen C, Li XR, Ran YY, Xu T, Zhang Y, Geng XK, Zhang Y, Du HS, Leak RK, Ji XM, Hu XM. Remote Ischemic Preconditioning - Mediated Neuroprotection against Stroke is Associated with Significant Alterations in Peripheral Immune Responses. *CNS Neurosci Ther*. 2016 Jan;22(1):43-52.
35. Zhang L, Xiong W, Xiong Y, Liu H, Liu Y. 17 β -Estradiol promotes vascular endothelial growth factor expression via the Wnt/ β -catenin pathway during the pathogenesis of endometriosis. *Mol Hum Reprod*. 2016 Jul;22(7):526-35.
36. Jiang L, Wang Y, Zhu D, Xue Z, Mao H. Alteration of histone H3 lysine 9 dimethylation in peripheral white blood cells of septic patients with trauma and cancer. *Mol Med Rep*. 2016 Dec;14(6):5467-5474.
37. Fan XY, Chen XY, Liu YJ, Zhong HM, Jiang FL, Liu Y. Oxidative stress-mediated intrinsic apoptosis in human promyelocytic leukemia HL-60 cells induced by organic arsenicals. *Sci Rep*. 2016 Jul 19;6:29865.
38. Xie Y, Wen H, Yan K, Wang S, Wang X, Chen J, Li Y, Xu Y, Zhong Z, Shen J, Chu D. *Toxoplasma gondii* GRA15II effector-induced M1 cells ameliorate liver fibrosis in mice infected with *Schistosomiasis japonica*. *Cell Mol Immunol*. 2016 May 9. doi: 10.1038/cmi.2016.21. [Epub ahead of print]
39. Li DY, Shi XJ, Li W, Sun XD, Wang GY. Ischemic preconditioning and remote ischemic preconditioning provide combined protective effect against ischemia/reperfusion injury. *Life Sci*. 2016 Apr 1;150:76-80.
40. Shi L, Zhang L, Li C, Hu X, Wang X, Huang Q, Zhou G. Dietary zinc deficiency impairs humoral and cellular immune responses to BCG and ESAT-6/CFP-10 vaccination in offspring and adult rats. *Tuberculosis (Edinb)*. 2016 Mar;97:86-96.
41. Liu ZJ, Chen C, Li XR, Ran YY, Xu T, Zhang Y, Geng XK, Zhang Y, Du HS, Leak RK, Ji XM, Hu XM. Remote Ischemic Preconditioning-Mediated Neuroprotection against Stroke is Associated with Significant Alterations in Peripheral Immune Responses. *CNS Neurosci Ther*. 2016 Jan;22(1):43-52.
42. Zhang T, Liu Y, Wang H, Zhang X, Zhu S, Xu P, Yin J, Ren G, Liu J, Li D. Recombinant NDV expressing cytokines or fliC confers a quick immune response against NDV challenge and resistance to maternal antibody. *Vet Microbiol*. 2016 Nov 30;196:107-117.
43. Jia H, Zhao T, Ji Y, Jia X, Ren W, Li C, Li M, Xiao Y, Wang H, Xu K. Combined nifuroxazide and SAT05f therapy reduces graft-versus-host disease after experimental allogeneic bone marrow transplantation. *Cell Death Dis*. 2016 Dec 1;7(12):e2507.
44. Du JY, Liang Y, Fang JF, Jiang YL, Shao XM, He XF, Fang JQ. Effect of systemic injection of heterogenous and homogenous opioids on peripheral cellular immuneresponse in rats with bone cancer pain: A comparative study. *Exp Ther Med*. 2016 Oct;12(4):2568-2576.
45. Wang Y, Qiu Y, Yin S, Zhang L, Shi K, Gao H, Zhang Z, He Q. A functional nanocarrier that copenetrates extracellular matrix and multiple layers of tumor cells for sequential and deep tumor autophagy inhibitor and chemotherapeutic delivery. *Autophagy*. 2017 Feb;13(2):359-370.
46. Zhang L, Xiong W, Li N, Liu H, He H, Du Y, Zhang Z, Liu Y. Estrogen stabilizes hypoxia-inducible factor 1 α through G protein-coupled estrogen receptor 1 in eutopic endometrium of endometriosis. *Fertil Steril*. 2017 Feb;107(2):439-447.

Version 2017.03.09